

Anti-Slc25A23 Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Bv35

Item Number: CM0015715

Introduction

Rabbit polyclonal antibody against SLC25A23, a mitochondrial ATP-Mg/Pi carrier regulating adenine nucleotide pools and calcium uptake. Validated for WB, IF/ICC, ELISA in human, mouse, and rat samples (52 kDa).

[Learn More](#)

Parameter	Value
Product Name	SLC25A23 Rabbit pAb
Remarks/Alias	APC2; MCSC2; SCAMC3; SCaMC-3; SLC25A23
Species	Human
Gene ID (Human)	79085
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-160 of human SLC25A23 (NP_077008.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BV35
Protein Weight	52kDa
Shipping	Ice bag