

Anti-Slc25A39 Polyclonal Antibody For Wb, If-P, Elisa - Q9Bzj4

Item Number: CM0020108

Introduction

Anti-SLC25A39 rabbit polyclonal antibody for WB, IF-P and ELISA applications, reactive with human, mouse and rat. Targets mitochondrial glutathione transporter, useful for redox biology and erythropoiesis research.

[Learn More](#)

Parameter	Value
Product Name	SLC25A39 Rabbit pAb
Remarks/Alias	CGI69; CGI-69; SLC25A39
Species	Human
GeneID (Human)	51629
GeneID	51629
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 170-250 of human SLC25A39 (NP_001137252.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZJ4
Protein Weight	39kDa
Shipping	Ice bag