


Anti-Slc40A1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Np59

Item Number: CM0014173



Introduction

Rabbit polyclonal antibody targeting human SLC40A1 (Ferroportin) validated for Western blot and ELISA. Shows cross-reactivity with human and mouse ferroportin. Suitable for iron homeostasis and hemochromatosis research.

[Learn More](#)

Parameter	Value
Product Name	SLC40A1 Rabbit pAb
Remarks/Alias	FPN; FPN1; HFE4; MTP1; IREG1; MST079; MSTP079; SLC11A3; SLC40A1
Species	Human
GeneID (human)	30061
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 200-270 of human SLC40A1 (NP_055400.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NP59
Protein Weight	63kDa
Shipping	Ice bag