


Anti-Slc6A6 Polyclonal Antibody For Wb And Elisa - P31641

Item Number: CM0013928



Introduction

Rabbit polyclonal antibody against human SLC6A6 (taurine transporter, TauT). Validated for Western blot (WB) and ELISA with cross-reactivity to human, mouse, and rat. Suitable for studying taurine transport, anti-aging mechanisms, and mitochondrial function.

[Learn More](#)

Parameter	Value
Product Name	SLC6A6 Rabbit pAb
Remarks/Alias	TAUT; HTRDC; SLC6A6
Species	Human
Gene ID (Human)	6533
Gene ID	6533
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 250-320 of human SLC6A6 (NP_001127839.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31641
Protein Weight	70kDa
Shipping	Ice bag