

Anti-Slc7A11/Xct Rabbit Polyclonal Antibody For Wb And Elisa - Q9Upy5

Item Number: CM0012833



Introduction

Rabbit polyclonal antibody against human SLC7A11/xCT, validated for WB and ELISA. Cross-reacts with mouse and rat. Unconjugated, 55 kDa. Ideal for oxidative stress, ferroptosis, and amino acid transport studies.

[Learn More](#)

Parameter	Value
Product Name	SLC7A11/xCT Rabbit pAb
Remarks/Alias	xCT; CCB1; SLC7A11/xCT
Species	Human
GeneID (Human)	23657
GeneID	23657
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SLC7A11/xCT (NP_055146.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UPY5
Protein Weight	55 kDa
Shipping	Ice bag