

Anti-Slfn11 Rabbit Polyclonal Antibody For Wb, Elisa - Q7Z7L1

Item Number: CM0020936



Introduction

SLFN11 Rabbit Polyclonal Antibody validated for WB, ELISA in human samples. Targets human SLFN11, a DNA replication inhibitor and antiviral restriction factor. Ideal for DNA damage and antiviral research.

[Learn More](#)

Parameter	Specification
Product Name	SLFN11 Rabbit pAb
Remarks/Alias	SLFN8/9; SLFN11
Species	Human
GeneID (Human)	91607
GeneID	91607
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 132-328 of human SLFN11 (NP_689483.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q7Z7L1
Protein Weight	103kDa
Shipping	Ice bag