

Anti-Slitrk6 Polyclonal Antibody For Wb, Elisa - Q9H5Y7

Item Number: CM0026480

Introduction

SLITRK6 Rabbit Polyclonal Antibody validated for WB and ELISA applications. Recognizes human, mouse, rat SLITRK6, a neurite outgrowth regulator linked to hearing and vision. Suitable for neurology and sensory research.

[Learn More](#)

Parameter	Value
Product Name	SLITRK6 Rabbit pAb
Remarks/Alias	DFNMYP; SLITRK6
Species	Human
GeneID (Human)	84189
GeneID	84189
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 621-816 of human SLITRK6 (NP_115605.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H5Y7
Protein Weight	95kDa
Shipping	Ice bag