


## Anti-Smac/Diablo Polyclonal Antibody For Wb, Ihc-P, Ip, Elisa - Q9Nr28

Item Number: CM0018948



### Introduction

Rabbit polyclonal antibody targeting human Smac/Diablo (SMAC; DFNA64) for WB, IHC-P, IP, ELISA. Cross-reacts with human, mouse, rat. Ideal for apoptosis and IAP interaction studies.

[Learn More](#)

| Parameter        | Value                                                                                                                    |
|------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Smac/Diablo Rabbit pAb                                                                                                   |
| Remarks/Alias    | SMAC; DFNA64; Smac/Diablo                                                                                                |
| Species          | Human                                                                                                                    |
| GeneID (human)   | 56616                                                                                                                    |
| GeneID           | 56616                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 56-239 of human Smac/Diablo (NP_063940.1). |
| Source           | Rabbit                                                                                                                   |
| Category         | Polyclonal Antibodies                                                                                                    |
| Application      | WB, IHC-P, IP, ELISA                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                        |
| SWISS            | Q9NR28                                                                                                                   |
| Protein Weight   | 27kDa                                                                                                                    |
| Shipping         | Ice bag                                                                                                                  |