


Anti-Smad1 Polyclonal Antibody For Wb, Elisa - Q15797

Item Number: CM0018307



Introduction

Rabbit polyclonal antibody targeting human Smad1 (aa 20-240). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Suited for BMP signaling and developmental biology research.

[Learn More](#)

Parameter	Value
Product Name	Smad1 Rabbit pAb
Remarks/Alias	BSP1; JV41; BSP-1; JV4-1; MADH1; MADR1; Smad1
Species	Human
GeneID (Human)	4086
GeneID	4086
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-240 of human Smad1 (NP_001003688.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15797
Protein Weight	52kDa
Shipping	Ice bag