


Anti-Smad4 Polyclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - Q13485

Item Number: CM0025045



Introduction

Rabbit polyclonal antibody against human Smad4 (SWISS: Q13485) for WB, IF/ICC, IHC-P, ELISA. Validated in human, mouse, rat. Key reagent for TGF-beta signaling and tumor suppression research.

[Learn More](#)

Parameter	Value
Product Name	Smad4 Rabbit pAb
Remarks/Alias	JIP; DPC4; MADH4; MYHRS; Smad4
Species	Human
GeneID (Human)	4089
GeneID	4089
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 320-552 of human Smad4 (NP_005350.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13485
Protein Weight	60 kDa
Shipping	Ice bag