


Anti-Smarcad1 Monoclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - Q9H4L7

Item Number: CM0005410



Introduction

Rabbit monoclonal antibody targeting human SMARCAD1, validated for WB, IP, IF/ICC, IHC-P, and ELISA. Detects SMARCAD1 in human, mouse, and rat samples. Ideal for chromatin remodeling and DNA repair studies.

[Learn More](#)

Parameter	Value
Product Name	SMARCAD1 Rabbit mAb
Remarks/Alias	HRZ; TYS; ETL1; HEL1; ADERM; BASNS
Species	Human
GeneID (Human)	56916
GeneID	56916
Immunogen	Recombinant protein corresponding to amino acids 928-1026 of human SMARCAD1 (NP_064544.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H4L7
Protein Weight	68 kDa / 117 kDa
Shipping	Ice bag