

Anti-Smarcb1/Snf5 Monoclonal Antibody For Wb, If/Icc, Elisa - Q12824

Item Number: CM0006100



Introduction

KO-validated rabbit monoclonal antibody against human SMARCB1/SNF5 (BAF47, INI1), suitable for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat samples. Ideal for chromatin remodeling and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] SMARCB1/SNF5 Rabbit mAb
Remarks/Alias	RDT; CSS3; INI1; SNF5; Snr1; BAF47; INI-1; MRD15; RTPS1; Sfh1p; hSNFS; SNF5L1; SWNTS1; PPP1R144; F5
Species	Human
Gene ID (Human)	6598
Gene ID	6598
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12824
Protein Weight	44kDa
Shipping	Ice bag