

Anti-Smarcb1/Snf5 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - Q12824

Item Number: CM0030148



Introduction

Rabbit polyclonal SMARCB1/SNF5 antibody (Q12824) validated for WB, IHC-P, IF/ICC, IP, ELISA. Cross-reacts with human, mouse, rat. Detects 44kDa BAF47 protein, a core component of the SWI/SNF chromatin remodeling complex involved in cell cycle control and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	SMARCB1/SNF5 Rabbit pAb
Remarks/Alias	RDT; CSS3; INI1; SNF5; Snr1; BAF47; INI-1; MRD15; RTPS1; Sfh1p; hSNF5; SNF5L1; SWNTS1; PPP1R144
Species	Human
Gene ID	6598 (Human)
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 50-150 of human SMARCB1/SNF5 (NP_003064.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q12824
Protein Weight	44kDa
Shipping	Ice bag