


Anti-Smmhc/Myh11 Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P35749

Item Number: CM0004986



Introduction

Rabbit monoclonal antibody targeting human SMMHC/MYH11 (myosin-11), suitable for WB, IF-P, IHC-P, ELISA in human, mouse, and rat. Recombinant immunogen; 227 kDa.

[Learn More](#)

Parameter	Value
Product Name	SMMHC/MYH11 Rabbit mAb
Remarks/Alias	AAT4; FAA4; SMHC; SMMHC; VSCM2; SMMS-1; SMMHC/MYH11
Species	Human
Gene ID (Human)	4629
Gene ID	4629
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1160-1313 of human SMMHC/MYH11 (NP_002465.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P35749
Protein Weight	227kDa
Shipping	Ice bag