


Anti-Smn1 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q16637

Item Number: CM0019791



Introduction

High-quality rabbit polyclonal antibody against human SMN1 (Survival Motor Neuron 1). Validated for WB, IHC-P, and ELISA, with cross-reactivity to mouse and rat. Ideal for splicing and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	SMN1 Rabbit pAb
Remarks/Alias	SMA; SMN; SMA1; SMA2; SMA3; SMA4; SMA@; SMNT; BCD541; GEMIN1; TDRD16A; T-BCD541; SMN1
Species	Human
GeneID (Human)	6606
GeneID	6606
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 11-169 of SMN1 (NP_000335.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16637
Protein Weight	32kDa
Shipping	Ice bag