

Anti-Snape2 Polyclonal Antibody For Wb And Elisa - Q13487

Item Number: CM0022730

Introduction

Rabbit polyclonal antibody against SNAPC2, a key subunit of the snRNA-activating protein complex. Validated for western blot and ELISA, it recognizes human and mouse SNAPC2 (~36 kDa). Ideal for transcription regulation studies and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	SNAPC2 Rabbit pAb
Remarks/Alias	SNAP45; PTFDELTA; SNAPC2
Species	Human
GeneID (Human)	6618
GeneID	6618
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-100 of human SNAPC2 (NP_003074.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13487
Protein Weight	36kDa
Shipping	Ice bag