

Anti-Snrnp70 Rabbit Monoclonal Antibody For Wb And Elisa - P08621

Item Number: CM0007136



Introduction

Rabbit monoclonal antibody targeting SNRNP70, validated for western blotting and ELISA, with reactivity to human and mouse samples. Ideal for spliceosome assembly and U1 snRNP research.

[Learn More](#)

Parameter	Value
Product Name	[KD Validated] SNRNP70 Rabbit mAb
Remarks/Alias	RPU1; Snp1; U1AP; U170K; U1RNP; RNPUIZ; SNRP70; U1-70K
Species	Human
GeneID (human)	6625
GeneID	6625
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 201-300 of SNRNP70 (NP_003080.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P08621
Protein Weight	52kDa
Shipping	Ice bag