

Anti-Snrpa1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P09661

Item Number: CM0013216



Introduction

High-quality SNRPA1 Rabbit Polyclonal Antibody for Western blot, IHC-P, and ELISA. Validated in human, mouse, and rat samples. Ideal for spliceosome and pre-mRNA splicing studies.

[Learn More](#)

Parameter	Value
Product Name	SNRPA1 Rabbit pAb
Remarks/Alias	Lea1, UZA', SNRPA1
Species	Human
GeneID (Human)	6627
GeneID	6627
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human SNRPA1 (NP_003081.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09661
Protein Weight	28kDa
Shipping	Ice bag