


Anti-Snrpb Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P14678

Item Number: CM0023188



Introduction

Rabbit polyclonal antibody against SNRNPB (Sm B/B'), a core spliceosomal protein. Validated for Western blot, IHC-P, and ELISA. Cross-reacts with human, mouse, and rat samples. For splicing and RNA processing studies.

[Learn More](#)

Parameter	Value
Product Name	SNRNPB Rabbit pAb
Remarks/Alias	COD; CCMS; SNRNPB1; SmB/B'; Sm-B/B'; snRNP-B; SmB/SmB'; SNRNPB
Species	Human
GeneID (human)	6628
GeneID	6628
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SNRNPB (NP_003082.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P14678
Protein Weight	25kDa
Shipping	Ice bag