


Anti-Socs2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - O14508

Item Number: CM0023557



Introduction

Human SOCS2 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Targets suppressor of cytokine signaling 2, a key regulator of growth hormone and erythropoietin receptor signaling. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	SOCS2 Rabbit pAb
Remarks/Alias	CIS2, SSI2, Cish2, SSI-2, SOCS-2, STATI2, SOCS2
Species	Human
GeneID	8835
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human SOCS2 (NP_001257400.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	O14508
Protein Weight	22kDa
Shipping	Ice bag