

Anti-Sox10 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Ip, Elisa - P56693

Item Number: CM0005293

Introduction

Anti-SOX10 monoclonal antibody validated for WB, IHC-P, IF/ICC, IP, ELISA in human, mouse, rat. Ideal for studying CNS myelination and transcription regulation. P56693.

[Learn More](#)

Parameter	Value
Product Name	SOX10 Rabbit mAb
Remarks/Alias	DOM; WS4; PCWH; WS2E; WS4C; SOX-10; SOX10
Species	Human
GeneID (Human)	6663
GeneID	6663
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 350-450 of human SOX10 (P56693).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P56693
Protein Weight	31kDa/50kDa
Shipping	Ice bag