

Anti-Sox18 Polyclonal Antibody For Wb, Elisa - P35713

Item Number: CM0022192

Introduction

High-quality rabbit polyclonal antibody against human SOX18, suitable for WB and ELISA. Cross-reacts with mouse and rat. For cardiovascular development and lymphangiogenesis research.

[Learn More](#)

Parameter	Value
Product Name	SOX18 Rabbit pAb
Remarks/Alias	HLTS; HLTRS; SOX18
Species	Human
GeneID (human)	54345
GeneID	54345
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SOX18 (NP_060889.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P35713
Protein Weight	41kDa
Shipping	Ice bag