

Anti-Spag8 Rabbit Polyclonal Antibody For Wb And Elisa - Q99932

Item Number: CM0025393



Introduction

Rabbit anti-SPAG8 polyclonal antibody validated for WB, ELISA. Cross-reacts with mouse and rat. Targets sperm-associated antigen 8, a ciliary protein involved in motility and spermatogenesis. For research use.

[Learn More](#)

Parameter	Details
Product Name	SPAG8 Rabbit pAb
Remarks/Alias	SMP1, BS-84, CT142, HSD-1, SPAG3, CILD28, hSMP-1, SPAG8
Species	Human
Gene ID (Human)	26206
Gene ID	26206
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-420 of human SPAG8 (NP_758516.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q99932
Protein Weight	51kDa
Shipping	Ice bag