


Anti-Spam1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P38567

Item Number: CM0006310



Introduction

Anti-SPAM1 rabbit monoclonal antibody validated for WB, IHC-P, ELISA. Specific for human SPAM1 (hyaluronidase PH-20), cross-reacts with rat. Ideal for sperm-egg adhesion and fertilization research.

[Learn More](#)

Parameter	Specification
Product Name	SPAM1 Rabbit mAb
Remarks/Alias	HYA1; PH20; HYAL1; HYAL3; HYAL5; PH-20; SPAG15; HEL-S-96n; SPAM1
Species	Human
Gene ID	6677
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 61-223 of human SPAM1 (NP_694859.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Rat
SWISS	P38567
Protein Weight	58kDa
Shipping	Ice bag