

Anti-Spdef Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O95238

Item Number: CM0016471



Introduction

Rabbit polyclonal antibody against human SPDEF (PDEF, SAM pointed domain-containing Ets transcription factor). Validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with rat. Suitable for prostate cancer and hormone-regulated epithelial tissue research.

[Learn More](#)

Parameter	Value
Product Name	SPDEF Rabbit pAb
Remarks/Alias	PDEF; bA375E1.3; SPDEF
Species	Human
GeneID (human)	25803
GeneID	25803
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-335 of human SPDEF (NP_036523.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	O95238
Protein Weight	38kDa
Shipping	Ice bag