


Anti-Spi1 Polyclonal Antibody For Wb, Ihc-P, Elisa - P17947

Item Number: CM0019975



Introduction

High-quality rabbit polyclonal antibody targeting human SPI1/PU.1 transcription factor. Suitable for WB, IHC-P, ELISA applications. Detects human and mouse SPI1. Essential for research on hematopoietic cell fate and gene regulation.

[Learn More](#)

Parameter	Value
Product Name	SPI1 Rabbit pAb
Remarks/Alias	OF; PU.1; AGM10; SFPI1; SPI-1; SPI-A; SPI1
Species	Human
GeneID (Human)	6688
GeneID	6688
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 17-98 of human SPI1 (NP_003111.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P17947
Protein Weight	31kDa
Shipping	Ice bag