

Anti-Spop Rabbit Monoclonal Antibody For Wb, Elisa - O43791

Item Number: CM0002452



Introduction

Recombinant rabbit monoclonal antibody targeting human SPOP (O43791). Validated for WB and ELISA; suitable for research on ubiquitination, proteasomal degradation, host-virus interactions, and DNA replication.

[Learn More](#)

Parameter	Value
Product Name	SPOP Rabbit mAb
Remarks/Alias	TEF2; BTBD32; NSDVS1; NSDVS2; NEDMACE; NEDMIDF; SPOP
Species	Human
Gene ID (Human)	8405
Gene ID	8405
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human SPOP (O43791).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43791
Protein Weight	42kDa
Shipping	Ice bag