

Anti-Spz1 Polyclonal Antibody For Wb And Elisa - Q9Bxg8

Item Number: CM0023611

Introduction

Rabbit polyclonal antibody against human SPZ1, validated in WB and ELISA. Recognizes human and rat SPZ1, a transcription factor involved in spermatogenesis. Ideal for reproductive biology and oncology research.

[Learn More](#)

Parameter	Value
Product Name	SPZ1 Rabbit pAb
Remarks/Alias	CILD38; NYD-TSP1; PPP1R148; SPZ1
Species	Human
GeneID (Human)	84654
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-430 of human SPZ1 (NP_115956.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q9BXG8
Protein Weight	49kDa
Shipping	Ice bag