

Anti-Sqstm1 Rabbit Polyclonal Antibody For Wb, Elisa - Q64337

Item Number: CM0022823



Introduction

Rabbit polyclonal antibody targeting SQSTM1/p62. Suitable for Western blot (WB) and ELISA. Detects human, mouse, and rat SQSTM1 with an expected band near 48 kDa. Ideal for autophagy and cell signaling studies.

[Learn More](#)

Parameter	Value
Product Name	SQSTM1 Rabbit pAb
Remarks/Alias	Osi, p62, A170, STAP, OSF-6, STONE14, SQSTM1
Species	Mouse
Gene ID	18412 (Mouse)
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of mouse SQSTM1 (NP_035148.1).
Source	Rabbit
Category	Polyclonal Antibody
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q64337
Protein Weight	48kDa
Shipping	Ice bag