

Anti-Sr-Bi Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q8Wtv0

Item Number: CM0002956

Introduction

Rabbit monoclonal antibody against human SR-BI (Q8WTV0). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Ideal for lipoprotein metabolism, viral entry, and phagocytosis research.

[Learn More](#)

Parameter	Specification
Product Name	SR-BI Rabbit mAb
Remarks/Alias	CLA1; SRB1; CLA-1; SR-BI; CD36L1; HDLCQ6; HDLQTL6
Species	Human
Gene ID	949
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human SR-BI (Q8WTV0)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTV0
Protein Weight	61kDa
Shipping	Ice bag