


Anti-SrpK2 Rabbit Polyclonal Antibody For Wb And Elisa - P78362

Item Number: CM0021376



Introduction

Rabbit polyclonal antibody targeting human SRPK2 (SRSF protein kinase 2). Validated for Western blot and ELISA in human, mouse, and rat. SRPK2 regulates mRNA splicing, apoptosis, and viral replication.

[Learn More](#)

Parameter	Value
Product Name	SRPK2 Rabbit pAb
Remarks/Alias	SFRSK2; SRPK2
Species	Human
GeneID (Human)	6733
GeneID	6733
Immunogen	Recombinant protein matching a sequence corresponding to amino acids 420-520 of human SRPK2 (NP_872633.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P78362
Protein Weight	78kDa
Shipping	Ice bag