

Anti-Srsf6 (Srp55) Rabbit Polyclonal Antibody For Wb And Elisa - Q13247

Item Number: CM0017936



Introduction

High-quality rabbit polyclonal antibody against human SRSF6 (SRp55) validated for WB and ELISA. Detects endogenous SRSF6 in human, mouse, and rat samples. Ideal for splicing research.

[Learn More](#)

Parameter	Value
Product Name	SRSF6 Rabbit pAb
Remarks/Alias	B52; SFRS6; SRP55; HEL-S-91; SRSF6
Species	Human
Gene ID (Human)	6431
Gene ID	6431
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human SRSF6 (NP_006266.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13247
Protein Weight	40kDa
Shipping	Ice bag