

Anti-Ssr2 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P43308

Item Number: CM0021001



Introduction

Rabbit polyclonal antibody against SSR2 (TRAP-beta), suitable for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for ER protein retention research.

[Learn More](#)

Parameter	Value
Product Name	SSR2 Rabbit pAb
Remarks/Alias	TLAP; HSD25; TRAPB; TRAP-BETA; SSR2
Species	Human
Gene ID (Human)	6746
Gene ID	6746
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P43308
Protein Weight	20kDa
Shipping	Ice bag