


Anti-Sstr1 Rabbit Polyclonal Antibody For Wb, Elisa - P30872

Item Number: CM0026124



Introduction

SSTR1 Rabbit Polyclonal Antibody detects human and rat somatostatin receptor type 1 (SSTR1/SRIF-2) in WB and ELISA applications. Recombinant immunogen targets the N-terminal domain (aa 1-56). Molecular weight: 43kDa.

[Learn More](#)

Parameter	Value
Product Name	SSTR1 Rabbit pAb
Remarks/Alias	SS1R, SST1, SS1-R, SRIF-2, SS-1-R, SSTR1
Species	Human
Gene ID	6751 (human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-56 of human SSTR1 (NP_001040.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P30872
Protein Weight	43kDa
Shipping	Ice bag