


Anti-Star Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P49675

Item Number: CM0013324



Introduction

StAR (STARD1) rabbit polyclonal antibody validated for WB, IP, ELISA in human, mouse. Targets mitochondrial steroidogenic acute regulatory protein (32 kDa, SWISS P49675). Ideal for steroidogenesis research.

[Learn More](#)

Parameter	Value
Product Name	StAR Rabbit pAb
Remarks/Alias	STARD1; StAR
Species	Human
Gene ID (Human)	6770
Gene ID	6770
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 64-285 of human StAR (NP_000340.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49675
Protein Weight	32kDa
Shipping	Ice bag