


Anti-Stat5B Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P51692

Item Number: CM0016256



Introduction

STAT5B rabbit polyclonal antibody for Western blot, IHC, and ELISA. Suitable for human and mouse samples. Detects signal transducer and activator of transcription 5B (90 kDa) involved in cytokine signaling. - P51692

[Learn More](#)

Parameter	Value
Product Name	STAT5B Rabbit pAb
Remarks/Alias	STAT5; GHISID2; STAT5B
Species	Human
Gene ID (Human)	6777
Gene ID	6777
Immunogen	Synthetic peptide corresponding to amino acids 594-794 of human STAT5B (NP_003143.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P51692
Protein Weight	90kDa
Shipping	Ice bag