

Anti-Sting/Tmem173 Monoclonal Antibody For Wb, Ip, If/Icc, If-P, Ihc-P, Elisa - Q86Wv6

Item Number: CM0003714



Introduction

Rabbit monoclonal antibody against human STING/TMEM173 (ERIS, MITA) for western blot, immunofluorescence, and ELISA. Cross-reacts with human, mouse, and rat samples. Ideal for innate immunity and STING pathway research.

[Learn More](#)

Parameter	Value
Product Name	STING/TMEM173 Rabbit mAb
Remarks/Alias	ERIS; MITA; MPYS; SAVI; NET23; STING; hMITA; hSTING; TMEM173; STING-beta; STING/TMEM173
Species	Human
Gene ID (Human)	340061
Gene ID	340061
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 150-379 of human STING/TMEM173 (NP_938023.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q86WV6
Protein Weight	42kDa
Shipping	Ice bag