


Anti-Strada Rabbit Polyclonal Antibody For Wb, Elisa - Q7Rtn6

Item Number: CM0023259



Introduction

Rabbit polyclonal antibody against human STRADA (STE20-related kinase adapter protein alpha), suitable for Western blotting and ELISA. Detects human and mouse STRADA, with a predicted molecular weight of ~48 kDa.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] STRADA Rabbit pAb
Remarks/Alias	LYK5; PMSE; Stlk; STRAD; NY-BR-96; STRAD alpha; DA
Species	Human
GeneID (Human)	92335
GeneID	92335
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q7RTN6
Protein Weight	48kDa
Shipping	Ice bag