

## Anti-Stxbp3 Polyclonal Antibody For Wb And Elisa - 000186

Item Number: CM0030190

### Introduction

Rabbit polyclonal antibody targeting human STXBP3 (syntaxin-binding protein 3, MUNC18C) for western blot (WB) and ELISA. Validated cross-reactivity with human, mouse, and rat samples. Ideal for GLUT4 trafficking research.

[Learn More](#)

| Parameter        | Value                                                                                                                |
|------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Name     | STXBP3 Rabbit pAb                                                                                                    |
| Remarks/Alias    | PSP; MUNC18C; UNC-18C; MUNC18-3; STXBP3                                                                              |
| Species          | Human                                                                                                                |
| GeneID (human)   | 6814                                                                                                                 |
| GeneID           | 6814                                                                                                                 |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 343-592 of human STXBP3 (NP_009200.2). |
| Source           | Rabbit                                                                                                               |
| Category         | Polyclonal Antibodies                                                                                                |
| Application      | WB, ELISA                                                                                                            |
| Cross Reactivity | Human, Mouse, Rat                                                                                                    |
| SWISS            | O00186                                                                                                               |
| Protein Weight   | 68kDa                                                                                                                |
| Shipping         | Ice bag                                                                                                              |