

Anti-Suclg1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P53597

Item Number: CM0018041



Introduction

Rabbit polyclonal antibody targeting SUCLG1 (succinyl-CoA ligase subunit alpha). Validated for WB, IF/ICC, ELISA in human, mouse, and rat. Suitable for TCA cycle and mitochondrial metabolism research.

[Learn More](#)

Parameter	Value
Product Name	SUCLG1 Rabbit pAb
Remarks/Alias	GALPHA; MTDPS9; SUCLA1; SUCLG1
Species	Human
Gene ID (Human)	8802
Gene ID	8802
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 41-346 of human SUCLG1 (NP_003840.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P53597
Protein Weight	36kDa
Shipping	Ice bag