

Anti-Sult1A2 Polyclonal Antibody For Wb, Elisa - P50226

Item Number: CM0017912



Introduction

High-quality rabbit polyclonal antibody targeting human SULT1A2 sulfotransferase, suitable for Western blot and ELISA applications. Detects endogenous SULT1A2 in human, mouse, and rat samples, supporting research in drug metabolism and cancer susceptibility.

[Learn More](#)

Parameter	Value
Product Name	SULT1A2 Rabbit pAb
Remarks/Alias	STP2; HAST4; P-PST; ST1A2; TSPST2; P-PST 2; SULT1A2
Species	Human
GeneID (Human)	6799
GeneID	6799
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 201-295 of human SULT1A2 (NP_001045.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50226
Protein Weight	34kDa
Shipping	Ice bag