


Anti-Supt16H (Spt16) Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Y5B9

Item Number: CM0019952



Introduction

High-quality SUPT16H (SPT16) rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat FACT complex subunit SPT16. Ideal for chromatin research.

[Learn More](#)

Parameter	Value
Product Name	SUPT16H (SPT16) Rabbit pAb
Remarks/Alias	CDC68; SPT16; NEDDFAC; FACTP140; SPT16/CDC68; SUPT16H (SPT16)
Species	Human
GeneID (Human)	11198
GeneID	11198
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-695 of human SUPT16H (SPT16) (NP_009123.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9Y5B9
Protein Weight	120kDa
Shipping	Ice bag