

Anti-Survivin Rabbit Monoclonal Antibody For Wb And Elisa - O15392

Item Number: CM0012608



Introduction

High-quality Anti-Survivin Rabbit Monoclonal Antibody validated for Western blot and ELISA. Exhibits cross-reactivity with mouse and rat. Ideal for apoptosis and cell proliferation research.

[Learn More](#)

Parameter	Value
Product Name	Survivin Rabbit mAb
Remarks/Alias	API4; EPR-1; Survivin
Species	Human
GeneID (human)	332
GeneID	332
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human Survivin (O15392).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O15392
Protein Weight	16kDa
Shipping	Ice bag