

Anti-Synpo (Synaptopodin) Polyclonal Antibody For Wb, Elisa - Q8N3V7

Item Number: CM0027301

Introduction

Rabbit polyclonal antibody targeting human SYNPO (Synaptopodin) validated for WB and ELISA. Detects endogenous levels in human, mouse, and rat. Ideal for dendritic spine and podocyte research.

[Learn More](#)

Parameter	Value
Product Name	SYNPO Rabbit pAb
Remarks/Alias	SYNPO
Species	Human
GeneID (Human)	11346
GeneID	11346
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human SYNPO (NP_001159680.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N3V7
Protein Weight	99kDa
Shipping	Ice bag