

Anti-Syntrophin Alpha 1 Monoclonal Antibody For Wb And Elisa - Q13424

Item Number: CM0003651

Introduction

Rabbit monoclonal antibody targeting human Syntrophin alpha 1 (Q13424). Validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for dystrophin complex and neuromuscular synapse research.

[Learn More](#)

Parameter	Value
Product Name	Syntrophin alpha 1 Rabbit mAb
Remarks/Alias	SNT1; LQT12; TACIP1; dj1187j4.5; Syntrophin alpha 1
Species	Human
GeneID (Human)	6640
GeneID	6640
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Syntrophin alpha 1 (Q13424)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13424
Protein Weight	54kDa
Shipping	Ice bag