

Anti-Tab3 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N5C8

Item Number: CM0018991

Introduction

High-quality rabbit polyclonal antibody targeting human TAB3 (NAP1/MAP3K7IP3). Validated for WB, IHC-P, and ELISA applications. Cross-reacts with mouse and rat. Ideal for studying NF-κB and JNK signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	TAB3 Rabbit pAb
Remarks/Alias	NAP1, MAP3K7IP3, TAB3
Species	Human
GenelD (Human)	257397
GenelD	257397
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 390-480 of human TAB3 (NP_690000.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N5C8
Protein Weight	79kDa
Shipping	Ice bag