


Anti-Tada2A Polyclonal Antibody For Wb, Elisa - 075478

Item Number: CM0030534



Introduction

Rabbit polyclonal antibody targeting human TADA2A (ADA2). Validated in WB and ELISA with cross-reactivity to mouse and rat. Suitable for chromatin remodeling and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	TADA2A Rabbit pAb
Remarks/Alias	ADA2, ADA2A, KL04P, hADA2, TADA2L, TADA2A
Species	Human
Gene ID	6871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-443 of human TADA2A (NP_001159577.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O75478
Protein Weight	52kDa
Shipping	Ice bag