

Anti-Tbl1Xr1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Bzk7

Item Number: CM0026697



Introduction

Anti-TBL1XR1 rabbit polyclonal antibody validated for WB and ELISA, targeting human and mouse TBL1XR1 (Gene ID: 79718, SWISS: Q9BZK7). Suitable for nuclear receptor transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	TBL1XR1 Rabbit pAb
Remarks/Alias	C21; DC42; IRA1; MRD41; TBLR1; TBL1XR1
Species	Human
Gene ID	79718
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human TBL1XR1 (NP_078941.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9BZK7
Protein Weight	56kDa
Shipping	Ice bag