

Anti-Tbx3 Rabbit Polyclonal Antibody For Wb, Elisa - 015119

Item Number: CM0025941



Introduction

Rabbit polyclonal antibody targeting human TBX3, validated for Western blot (WB) and ELISA. Cross-reacts with mouse. Designed for studies in developmental biology, transcriptional regulation, and T-box factor research - 015119.

[Learn More](#)

Parameter	Value
Product Name	TBX3 Rabbit pAb
Remarks/Alias	UMS; XHL; TBX3-ISO; TBX3
Species	Human
GeneID (Human)	6926
GeneID	6926
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-480 of human TBX3 (NP_057653.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	015119
Protein Weight	79kDa
Shipping	Ice bag