


Anti-Tdo2 Polyclonal Antibody For Wb, If/Icc, Elisa - P48775

Item Number: CM0015297



Introduction

Rabbit polyclonal antibody against human TDO2 (tryptophan 2,3-dioxygenase), validated for WB, IF/ICC, and ELISA. Detects endogenous TDO2 in human, mouse, and rat samples. Ideal for tryptophan metabolism research.

[Learn More](#)

Parameter	Value
Product Name	TDO2 Rabbit pAb
Remarks/Alias	TO; TDO; TPH2; TRPO; HYPTRP; TDO2
Species	Human
GeneID (Human)	6999
GeneID	6999
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 100-200 of human TDO2 (NP_005642.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P48775
Protein Weight	48kDa
Shipping	Ice bag